Tensoval[®]

compact

Instructions for use

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www.hartmann.info

(0210)

 Please read this introduction carefully before use as correct measurement can only be carried out following these instructions. Should you have any further questions on our product or on measuring blood pressure, please contact our Customer Service (a list of addresses are given in the warranty document). On the Internet, you can visit our homepage on www.hartmann.info.

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Congratulations on your purchase of this HARTMANN quality product. Tensoval compact is a semi-automatic blood pressure monitor which uses intelligent measuring technology to take rapid and reliable measurement of systolic and diastolic blood pressure as well as pulse rate using an oscillometric measuring method.

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
1. General information on blood pressure

The human heart beats approx. 60 to 80 times per minute. This causes blood to be pumped into the arterial vascular system and supplies the body with oxygen and necessary nutrients. A constant pressure known as blood pressure is required to enable the blood to flow into the smallest of blood vessels. In measuring blood pressure, two values are determined. The systolic value is the maximum pressure as the heart beats. The diastolic value is the minimum pressure between two

heartbeats. Blood pressure is expressed as, e.g. 120/80 in millimetres of a mercury column (mmHg).

Blood pressure changes constantly in each individual and thus creates the requirements necessary for the functional capacity of the body. Thus fluctuations in blood pressure are perfectly normal. However, if blood pressure values are constantly elevated at rest, hypertension or high blood pressure is present.

The World Health Organisation (WHO) has set out the following limit values:

Assessment	Systolic pressure	Diastolic pressure
Optimal	up to 120 mmHg	up to 80 mmHg
Normal	up to 130 mmHg	up to 85 mmHg
Normal limit values	130 – 139 mmHg	85 – 89 mmHg
Grade 1 hypertension	140 – 159 mmHg	90 – 99 mmHg
Grade 2 hypertension	160 – 179 mmHg	100 – 109 mmHg
Grade 3 hypertension	over 180 mmHg	over 110 mmHg

High blood pressure is one of the most common causes of disability and death. Yet many people do not know they are suffering from high blood pressure. This is because the problem is only noticed when the conditions it leads to are at an advanced stage. These may be life-threatening, such as heart attack, stroke or kidney failure. Only regular blood pressure monitoring enables high blood pressure to be detected early. You have chosen blood pressure self-measurement and this will contribute to your health care. This is the only way to ensure efficient protection against the dangerous results which may affect your heart and circulation.

2. Important notes for self-measurement

- Self measurement of blood pressure does not replace either treatment or regular check-ups by your doctor. The drug dose prescribed by your doctor may therefore on no account be altered without consultation.
- Self-measurement of blood pressure must be conducted regularly

and over a long period of time. With the assistance of the blood pressure values recorded regularly by you in your blood pressure diary, high blood pressure can be detected early and your doctor can initiate optimal treatment. Individual readings are situation-related and thus are not useful.

- In order to compare values, it is important that you always take your blood pressure at rest. So you should relax for 5 minutes before taking a measurement.
- Even slight changes in internal and external factors (e.g. deep breathing, stimulants, talking, excitement, climatic factors) lead to fluctuations in blood pressure. This is why your doctor and pharmacist often obtain different readings.
- Measure your blood pressure regularly in the morning and evening at about the same time, as blood pressure fluctuates during the course of the day.
- You should not move or talk during the measurement process.
- Always measure blood pressure on the naked upper arm, and

make sure that items of clothing rolled up do not impair circulation. Always measure blood pressure on the same arm and rest the arm, relaxed on a support.

- The cuff is suitable for an upper arm circumference of between 22 and 32 cm. Outside these limits correct readings may no longer be guaranteed.
- When taking repeated measurements, relax for at least one minute between measurements as otherwise false results are obtained.
- In severe cardiac arrhythmias (cases of disturbed heart rhythm), measurements should only be taken in consultation with a doctor. Due to the oscillometric measuring method, measurements taken with this device may in some instances be false or fail to produce a reading (Err).
- Monitoring blood pressure is extremely important during pregnancy, as blood pressure may be altered by the pregnancy. The results should, however, only be interpreted in consultation with the doctor.

3. Control displays and symbols

-  Flashes when the device is measuring and the pulse is being taken
-  Change battery
-  The cuff is not adequately inflated. Pump again to raise the set pressure value.
-  Deflate the cuff fully by pressing the exhaust button.
-  Measuring error, cf. Chap. 11
-  Protection from electric shock (type B)
-  Follow the operating instructions
-  Pumping bulb contains natural rubber latex which can cause allergic reactions.


4. Current supply

Open the battery cover on the underside of the device. Insert batteries (see chap. 12, Technical specifications), ensuring correct polarity (“+” and “-”). Close the battery lid.

- If you do not intend to use the device for a longer period of time, please remove the batteries.
- In the interests of environmental protection, used batteries should not be disposed of in household waste. Please observe current regulations when disposing of the batteries or use municipal collecting bins.
- Do not discard with household waste, but use specially designated disposal facilities (recycling centres etc.).


5. Connecting the tubes

6. Applying the cuff

- Measurement should be taken on the naked arm which usually has the higher blood pressure value.
- When using the draw-clamp cuff, thread the end of the cuff through the metal hoop if necessary, so as to form a loop. In this instance the Velcro fastener must be on the outside. Wrap the cuff round the upper arm, the tube lies in the centre of the elbow bend, pointing towards the hand. Take the free end of the cuff, wrap it firmly round the arm and close the Velcro fastener.
- Check that the lower edge of the

cuff is approx. 2.5 cm away from the bend of the elbow (see fig.).

- The cuff should be firm but not too tight. You should be able to push one or two fingers between the arm and the cuff. Note: Uneven wrapping of the cuff can lead to inaccurate results.

7. Measuring blood pressure


- It is recommended that blood pressure measurements are carried out whilst seated. Place your arm with the palm upwards on a support and make sure the cuff is at the level of the heart.
- Do not turn the device on until the cuff has been applied, otherwise the cuff can become damaged through the resulting excess pressure.
- Press the START/STOP button. All display symbols appear followed by the last measured value. A

flashing arrow pointing downwards shows that the device is checking itself automatically. A beep sounds and the “0” appears in the display which signals that the unit is now ready to take a new measurement. If the downward arrow symbol “▼” is displayed, there is some air remaining in the cuff. Release air using the exhaust button until the “0” appears.


- Pressurize the cuff using the pumping bulb until you have a pressure value which is suitable for you. It should be approximately 30 mmHg higher than your regular systolic blood pressure and is shown on the display. When pressurization is completed, Tensoval compact automatically begins the measurement. If the cuff was not inflated enough or if the measurement is disturbed, the arrow pointing upwards “▲” appears. Pump again until a high enough pressure is reached.
- Important: do not move or talk throughout the entire measuring procedure.
- As pressure in the cuff decreases,

the heart symbol and the falling cuff pressure are displayed. A beep sounds which is synchronised with your heartbeat.

- A beep indicates the end of measurement. The systolic and diastolic readings appear simultaneously on the display, with the pulse rate appearing alternately (see fig.)

- The reading is stored automatically.
- Deflate the cuff completely by pressing the exhaust button.

- To switch off the device, simply press the START/STOP button. Otherwise the device will switch itself off automatically after 3 minutes.
- If you wish to stop measurement for any reason, simply press the exhaust button. The measuring process is interrupted and an automatic fall in pressure occurs.

8. Maintenance

- Do not expose the device to extreme temperatures, humidity, dust or direct sunlight.
- This device consists of high quality precision components. Do not drop the device and protect it against hard knocks.
- Never open the device. Repairs may only be carried out by qualified individuals.
- Only use a soft, moistened cloth

to clean the device. Do not use detergents or solvents.

- The cuff can be cleaned carefully with a moist cloth and mild, pH-neutral soap. Do not completely immerse the cuff in water.

9. Checking the measurement function

We recommend that measurement function be checked at intervals of two years in the case of professional users, such as pharmacies, medical practices and clinics. You should also observe the national

regulations, such as, in Germany, the Decree on Sales of Medical Devices (Medizinprodukte-Betreiberverordnung). Measurement function checks can be carried out either by competent authorities, or authorised maintenance providers for a fee.

10. Warranty

This product comes with a 3-year warranty. Details of warranty conditions and contact addresses can be found in the separate warranty document.

11. Error displays

Error which has occurred	Possible causes
Device will not switch on	No batteries, they have been inserted incorrectly or are empty.
▲	The cuff was not inflated enough. Pump again until a suitably high pressure is reached.
Err 280	The cuff was inflated to more than 280 mmHg. Deflate the cuff by pressing the exhaust button.

Error which has occurred	Possible causes
Err	Cuff was applied incorrectly or too loosely. Measurement was disturbed due to speaking or moving.
PUL Err	The pulse could not be measured accurately.
	The battery is empty. Replace the battery.
--- ---	The device has a technical problem. Contact the shop where you bought the device or the HARTMANN representation for your country.
Cuff will not inflate	Cuff connection is incorrectly positioned in device.
Implausible reading	Incorrect cuff size. Cuff placed on top of clothing. Rolled up clothing impedes blood circulation. Moving, talking or excitement during measurement. Deep breathing during measurement process. No relaxation before measurement. Stimulants taken just before measurement.

Switch the device off if an error symbol appears. Check the possible causes and note the instructions on self-measurement in Chapter 2. Relax for a minute and take the measurement again. Do not move or talk during the measurement.

12. Technical specifications

Measuring method:	oscillometric
Display range:	0 – 280 mmHg
Measuring range:	Systole (SYS): 50 – 270 mmHg Diastole (DIA): 30 – 250 mmHg Pulse: 40 – 200 Pulses/min.
Technical accuracy of measurement:	cuff pressure: ± 3 mmHg Pulse: ± 5 % of displayed pulse frequency
Clinical accuracy of measurement:	complies with the EN 1060 standard, part 3.
Power supply:	4 x 1.5 V Mignon R6/LR6
Automatic cut-out:	3 minutes after end of measuring process.
Cuff:	Normal cuff 22 – 32 cm
Storage capacity:	1 measurement
Operating conditions:	Ambient temperature: $+10$ °C to $+40$ °C Relative humidity: 30 – 85 %
Storage/transport conditions:	Ambient temperature: -20 °C to $+60$ °C Relative humidity: 10 – 95 %
Serial number:	see battery compartment

13. Legal requirements and guidelines

Tensoval compact complies with the requirements of the EC directive, 93/42/EEC on Medical devices and bears the CE-marking.

The device complies with the EN 1060 standard relating to non-invasive blood pressure monitors Part 1: General requirements EN 1060-1: 1995 and Part 3: Additional requirements for electromechanical blood pressure measuring systems EN 1060-3: 1997. Clinical testing of measurement accuracy was performed according to the testing protocol recommended by EN 1060-3 (according to DIN 58130 and ANSI/AAMI SP10-1992).

Manufacturer:
PAUL HARTMANN AG
89522 Heidenheim
Germany

14. Accessories and spare parts

To ensure measurement accuracy, only use original HARTMANN accessories which may be obtained from your pharmacist or specialist medical supplier.

- Draw-clamp cuff
Code No.: 900 179
Upper arm circumference
22 – 32 cm
- Pumping bulb
Code No.: 900 190

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